

Which Podo autoclave is right for you?

Discover our wide range of autoclaves for podiatry practices.



Product highlights	 SteriHero Speed+ (Class B)	 SteriHero Podo 18 S (Class S)	 SteriHero Podo 18 (Class B)	 Vacuclave 318 (Class B)	 Vacuklav 41 B+ Evolution (Class B)
Cooling technology: for reliable operation in everyday practice	Stand-alone with air cooling	Stand-alone with air cooling	Stand-alone with air cooling	Stand-alone with air cooling	Stand-alone with air cooling
Water supply: with demineralized water for sterilization	External water tanks	Integrated water tanks or optional water treatment	Integrated water tanks or optional water treatment	Integrated water tanks or optional water treatment	Integrated water tanks or optional water treatment
Chamber sizes: for your individual loading requirements	5 liters	17 liters	17 liters	17 liters	18 liters
Maximum capacity: for loading with individual instruments	Up to 2 kg	Up to 4 kg	Up to 5 kg	Up to 7 kg	Up to 9 kg
Number of trays: for efficient use of the sterilization chamber	2 pcs.	5 pcs.	5 pcs.	5 pcs.	8 pcs.
Number of sterilization containers: for loading with MELAstore Boxes	✗	3x container (Attention: Observe maximum weight!)	3x container	3x container	4x container
Cycle times: for one sterilization process (without drying)*	6.5 – 16 min	13 – 23 min	15 – 27 min	15 – 27 min	8 – 21 min
Smart-Touch Display: for intuitive operation	4.3"	4.3"	4.3"	4.3"	5.7"
Integrated documentation interfaces: for protocol storage via network or data carrier	✓	✓	✓	✓	✓
ProControl: Batch approval and user authentication directly on the device	✗	✗	✗	✓	✓
Printing of barcode labels: with the help of an optional label printer	✗	✗	✗	✓	✓
DRYtelligence: for optimum drying results and a reduction in drying time of up to 80 %	✗	✗	✗	✓	✓
Energy-saving function: to deactivate the heating and the display between sterilization processes	✓	✓	✓	✓	✓
Start time preselection: for an automatic program start at the desired time	✗	✗	✗	✗	✓
Maintenance cycles: to ensure long-term quality and reliability	4.000 cycles or 2 years	2.000 cycles or 2 years	2.000 cycles or 2 years	2.000 cycles or 2 years	2.000 cycles or 2 years

* Time span for minimum load in the Quick S program to maximum load in the Universal B program

The Podo Autoclave Guide:

Podiatry practices face the same challenges as doctors' and dentists' practices: compliance with strict hygiene regulations according to the Robert Koch Institute. However, comparable financial resources are often lacking. SteriHero takes on this task and enables efficient hygiene processes with a top price-performance ratio.



Maximum reliability:

Thanks to the unique quality – 100% Made in Germany



Top price-performance ratio for podiatry:

With smart-touch display and integrated documentation



Excellent support:

Through your direct MELAG contact persons and our customer service for maintenance, validation and service



SteriHero Speed+

Groundbreaking performance, an unbelievably straightforward operating concept and extreme flexibility for podiatry practices: the fastest compact autoclave in its class allows sterilization of your instruments between treatments from 6.5 minutes (plus drying).



SteriHero Podo 18 S

No matter what big tasks the day has in store for you, SteriHero Podo 18 has your back. That's because the brand new class S stand-alone autoclave now also meets the requirements for a medical device and is perfectly tailored to podiatry: from the functions to the price!



SteriHero Podo 18

Your new superpower: the price-performance hero with 17 liters of volume offers the safest sterilization thanks to the fractionated pre- and post-vacuum process of class B - even for complex hollow instruments in your podiatry practice.

